

Overview of Rapid Antigen Test Specifications

	Abbott Panbio	BD Veritor	BTNX
Specimen Collection	Nasal	Nasal	Nasal or NP
Sensitivity	91.1%	83.9 - 93.5%	90.2%
Specificity	99.7%	99.3%	100%
Storage Temperature		2-30 °C	
Storage Requirements	Smaller box	Larger box	Smaller box
Testing TAT		15 mins	
Throughput (Analyzer)	No analyzer	No Analyzer <i>'analyzer not required, but can still be used</i>	No Analyzer
Shelf Life	12 months Regulator to advise on extending expiry date shortly	12 months Regulator to advise on 3-4 month shelf-life extension shortly	24 months
Ease of At-Home Testing	1 bottle of buffer; must be titrated into tubes	Pre-filled tubes, in sealed bags of 10	1 bottle of buffer; must be titrated into tubes